



**BOSCO**  
MEDICAL AUSTRALIA

# Salter Labs



## SOFT Nasal Oxygen Cannula

### Designed to deliver comfort at every point that touches the patient

No other cannula ever felt like this. The new 16SOFT nasal oxygen cannula from Salter was carefully designed to provide a higher level of comfort for your nasal oxygen patients. Comfortable patients can be happy patients, who follow their treatments, just by keep their cannulas on. The softness is everywhere. Curved, comfortable prongs fit gently into patients nostrils without pressing too hard upon the upper lip. Soft, pliable, stretchable tubing curves around the cheeks and ears, lessening the pressure where most device-related pressure ulcers develop.



It features a totally redesigned slide bolo to resist twisting and turning. Every inch of the new 16SOFT nasal oxygen cannula from Salter is designed to bring a touch of comfort and softness wherever it touches a patients face. No wonder Salters 1600 series oxygen cannula is most often worn by long term oxygen patients.


**LATEX  
FREE**




In The Hospital



In The Home



In Nursing Care



### Patient Benefits

- 1 Soft, 3-channel safety tubing minimizes kinking and disturbed oxygen flow
- 2 Pliable, stretchable tubing feels soft and gentle to the touch
- 3 Unique bolo design helps resist constant twisting
- 4 Curved, soft prongs avoid feeling bothersome and intrusive
- 5 Smooth, curved face-piece conforms to natural anatomical facial structure

| Cannula Type    | Weight  | LPM         | Line Length | QTY    | Product Code    |
|-----------------|---------|-------------|-------------|--------|-----------------|
| Premature       | ≤1400g  | Up to 3LPM  | 7' (2.1m)   | Box 50 | 16SOFT-PRE-7-50 |
| Neonate         | >1400g  | Up to 3LPM  | 7' (2.1m)   | Box 50 | 16SOFT-NEO-7-50 |
| Infant          | >2.5kg  | Up to 3LPM  | 7' (2.1m)   | Box 50 | 16SOFT-INF-7-50 |
| Toddler         | >4.0kg  | Up to 3LPM  | 7' (2.1m)   | Box 50 | 16SOFT-TOD-7-50 |
| Paediatric      | >12.5kg | Up to 6LPM  | 7' (2.1m)   | Box 50 | 16SOFT-PED-7-50 |
| Adult           | -       | -           | -           | Box 50 | 16SOFT-0-50     |
| Adult           | -       | Up to 6LPM  | 7' (2.1m)   | Box 50 | 16SOFT-7-50     |
| Adult           | -       | Up to 6LPM  | 14' (4.2m)  | Box 50 | 16SOFT-14-50    |
| Adult           | -       | Up to 6LPM  | 25' (7.6m)  | Box 50 | 16SOFT-25-25    |
| Adult High Flow | -       | Up to 15LPM | 7' (2.1m)   | Box 50 | 16SOFT-HF-7-25  |
| Adult High Flow | -       | Up to 15LPM | 14' (4.2m)  | Box 50 | 16SOFT-HF-14-25 |