



1600HF - High Flow Nasal Cannula

Salter design innovation. The high-performance, high-flow nasal cannula.

Wear-ability and comfort are standard with Salter cannulas. The comfort and convenience of oxygen therapy with nasal cannulas is now readily available for patients requiring higher oxygen flows. Give them the freedom of Salter's High Flow Nasal Cannula.

The Salter High Flow Cannula can replace cumbersome face masks with its enhanced reservoir face-piece and 3-channel safety tubing for flow rates up to 15 L/min.

Extensive testing provides ample evidence of comparative performance.



Salter features

- Delivers O₂ flows up to 15 L/min with 3-channel tubing
- Unique headset tubing with patented low-profile slide bolo
- Simple to maintain and operate dry bubble humidifier designed specifically for higher O₂ delivery

Problems solved

- Eliminates the need to use face masks
- Increased mobility
- Less cumbersome and claustrophobic than a mask
- Quieter operation than heated high-flow cannula systems
- Simple, secure adjustments

Product Code	Description	QTY
1600HF	Adult, Clear cannula with enhanced reservoir face-piece with 7' (2.13m) Bright Green supply tubing	Box 25