



BOSCO
MEDICAL AUSTRALIA

Aspen Surgical



Dr. Fog Anti-Fog Bottle with Sponge

WHEN IT COMES TO SURGICAL PROCEDURES, INSIST ON CLARITY.

Anti-Fog products for uninterrupted endoscopic procedures.

The Dr. Fog Anti-Fog brand is unique because it is one of the original brands of Anti-Fog on the market. With this trusted Dr. Fog brand, Aspen Surgical is a market share leader in the Anti-Fog market. DF-3100 the original anti-fog solution comes with an easy to squeeze bottle and sponge kit.



**LATEX
FREE**

- Intended as an anti-fog agent for endoscopic camera lenses.
- Removes fog formation on lens and allows the procedure to be uninterrupted.
- Adhesive backed sponge allows for secure, safe, and convenient placement in sterile field.
- Sponge is radiopaque.
- Sterile.
- Not made with natural rubber latex.

Product Code	Description	QTY
DF-3100	Original Solution With Sponge	Box 20



BOSCO
MEDICAL AUSTRALIA

Bosco Medical Australia Pty Ltd
Unit 23, 160 Lytton Road
Morningside, QLD 4170

P: + 61 7 3399 1799 F: + 61 7 3395 6736
Email: orders@boscomed.com.au
Website: www.boscomed.com.au