




# Needle Counters

## Deluxe Needle Counters

- Full-surface magnet provides extra space for odd-shaped sharps, pins, blades or needles
- Numbers on magnet allow easy sharps placement and positioning
- Easy-to-read numbers and description facilitate correct count of blades and needles
- Red colour signifies biohazard waste
- Not made with natural rubber latex



### MAGNET

Product Code	Description	QTY
25-0104	Needle counter double magnet red (30 count 60 capacity)	Box 50

### MAGNET/FOAM COMBO


Product Code	Description	QTY
25-0402	Needle counter foam block/magnet red (10 count 20 capacity)	Box 25

### FOAM

Product Code	Description	QTY
25-0502	Needle counter foam block/adhesive red (20 count 40 capacity)	Box 20
25-0504	Needle counter foam block/adhesive red (40 count 80 capacity)	Box 12
25-0801	Needle counter foam strips red (30 count 60 capacity)	Box 50

## DEBLADE NEEDLE COUNTERS

- Red colour signifies biohazard waste
- Featuring “Lock Box” design
- Scalpel handle guide facilitates use of the blade remover
- Overlapping box edges prevent sharps from escaping when closed
- User passive locking latch engages automatically to prevent inadvertent opening
- Optional blade remover decreases risk of clinician injury



### MAGNET

Product Code	Description	QTY
25-1104Y	DeBlade Needle Counter, Double Magnet yellow (60 Count 120 Capacity)	Box 50

### MAGNET/FOAM COMBO

Product Code	Description	QTY
25-1401	DeBlade Needle Counter, Foam Block/Magnetic red (20 Count 40 Capacity)	Box 20
25-1401Y	DeBlade Needle Counter, Foam Block/Magnetic yellow (20 Count 40 Capacity)	Box 20

### FOAM


Product Code	Description	QTY
25-1701	DeBlade Needle Counter, Foam Strips/Magnetic red (15 Count 30 Capacity)	Box 25



# Needle Counters

## CLEAR TOP NEEDLE COUNTERS

- Full-surface magnet provides extra space for odd-shaped sharps, pins, blades or needles
- Numbers on magnet allow easy sharps placement and positioning
- Easy-to-read numbers and description facilitate correct count of blades and needles
- See-through cover available, allowing clear box visualization when closed
- Red colour signifies biohazard waste
- Featuring "Lock Box" design
- User passive locking latch engages automatically to prevent inadvertent opening
- Overlapping box edges prevent sharps from escaping when closed
- Not made with natural rubber latex



## CLEAR TOP

Product Code	Description	QTY
25-0103	Needle counter single magnet/clear top red (10 count 20 capacity)	Box 50
25-0103Y	Needle counter single magnet/clear top yellow (10 count 20 capacity)	Box 50
25-0802	Needle counter foam strips red (15 count 30 capacity)	Box 50