



Original Cannulas

Works with the anatomy rather than against it...A Salter Innovation.

The Original Salter-Style Cannula is for oxygen flows up to 6 LPM. The clear, lightweight adult cannula is made with the Salter-Style manufacturing process, producing a soft, anatomically designed, one piece, seamless cannula that is comfortable to wear.



Salter features

- Anatomically curved to fit upper lip
- Curved, tapered nasal prongs
- Over-the-ear style headset
- Three-channel safety tubing
- Soft, lightweight, clear material

Benefits

- Eliminates flap or ridge which can cause irritation
- Improved anatomical fit for increased comfort
- Secure positioning
- Uninterrupted flow of oxygen
- Better comfort for long-term use

Salter-Style Original Cannula - 1600

Clear, lightweight cannula with unique, one piece anatomical design eliminates the need for flap, which can cause irritation. The result is maximum patient comfort. Over-the-ear style head set ensures secure positioning, allowing freedom of movement.

Product Code	Description	QTY
1600	Salter Adult Cannula with 7' (2.13m) supply tubing, curved prongs, over the ears, low flow (up to 6 LPM)	Box 50
1600-1	Cannula, Adult with 1'(0.3m) supply tubing, curved prongs, over the ears. Oxygen flow up to 6LPM	Box 50
1606	Salter Adult Cannula without supply tubing, curved prongs, over the ears, low flow (up to 6 LPM)	Box 50