



ACCUBOARD - Disposable Neonatal, Infant and Pediatric Limb-Boards



Ameritus AccuBoards feature a formable core that gently conforms to a baby, allowing small limbs to be comfortably positioned and stabilized during I.V. procedures.

- An aluminum core that may be easily contoured to virtually any shape and position for maximum patient comfort and support
- Soft, lightweight fabric and durable foam pad provide a gentle touch to fragile skin
- Ameritus AccuBoards are also available with soft, foam straps for securing the board to the patient
- Ameritus Limb-Boards feature a soft, spun fabric covering, thick foam padding, and an inner core designed for additional comfort and support
- Small, easy and cost effective to use
- Latex Free and DEHP Free
- Made in the USA

LATEX FREE


DEHP-free

Product Number	Description	QTY
403ABV	Neonatal (With Straps) - 1.91 cm x 6.67 cm	Box 48
404ABV	Infant (With Straps) - 4.13 cm x 9.85 cm	Box 48
405ABV	Pediatric (With Straps) - 4.75 cm x 12.38 cm	Box 48