



Anaesthetic Masks - Single Use

- Face Mask devised for anaesthesia with inflatable cushion
- Masks available with choice of light strawberry scents
- Masks available in Premium or Standard
- Hook rings included & sizes correspond to colour-coding for easy identification
- Single patient use
- Mask made from PVC, hook ring made from Polycarbonate
- Latex Free

**LATEX
FREE**



Product Code	Product Description	Size	Hook Ring	QTY
B1000	Premium Pre Neonate Strawberry Scented	0	-	Box 50
B1001	Premium Neonate Strawberry Scented	1	-	Box 50
B1001US	Premium Neonate Unscented	1	-	Box 50
B1002	Premium Infant Strawberry Scented	2	-	Box 50
B1002US	Premium Infant Unscented	2	-	Box 50
B1003	Premium Paediatric Strawberry Scented	3	Yellow	Box 50
B1003US	Premium Paediatric Unscented	3	Yellow	Box 50
B1003.5	Premium Medium Paediatric Strawberry Scented	3.5	Clear	Box 50
B1003.75	Premium Large Paediatric Strawberry Scented	3.75	Purple	Box 50
B1003.75US	Premium Large Paediatric Unscented	3.75	Purple	Box 50
B1004	Premium Small Adult Unscented	4	Green	Box 50
B1004SB	Premium Small Adult Strawberry Scented	4	Green	Box 50
B1005	Premium Medium Adult Unscented	5	Red	Box 50
B1006	Premium Large Adult Unscented	6	Blue	Box 50
B1101	Standard Neonate Strawberry Scented	1		Box 50
B1102	Standard Infant Strawberry Scented	2		Box 50
B1103	Standard Paediatric Strawberry Scented	3	Yellow	Box 50
B1103.75	Standard Large Paediatric Strawberry Scented	3.75	Purple	Box 50
B1104	Standard Small Adult Unscented	4	Green	Box 50
B1105	Standard Medium Adult Unscented	5	Red	Box 50
B1106	Standard Large Adult Unscented	6	Light Blue	Box 50

