




Pennine Healthcare

Uterine Aspiration Curette System

FLEXI-CUFF ASPIRATION TUBING

- Sterile, double wrapped individual peel pouches
- Leak proof seals
- Easily connects to curettes and Karman catheters
- Loose male connector in each pack for connection to suction machines where required
- Made of high density polythene
- Flexi-Cuff provides extra flexibility at a crucial point offering ease of use and improved technique



Wide Bore Aspiration Tubing (Flexi-Cuff)


Product Number	Description	QTY
AEKHS3000	2m wide bore + rotary handle	Each
AEKHS3030	3m wide bore + rotary handle	Each

ASPIRATION CURETTES

- Clear, rigid plastic curettes for uterine aspiration
- Moulded markings for sizing and correct tip orientation
- Seven sizes available in either curved or straight designs
- Connects to Pennine rotary swivel handle and wide bore, high vacuum tubing

Product Code	Curved Tip	QTY
AEKC1007	7mm	Box of 50
AEKC1008	8mm	Box of 50
AEKC1009	9mm	Box of 50
AEKC1010	10mm	Box of 50
AEKC1011	11mm	Box of 50
AEKC1012	12mm	Box of 50

Product Code	Straight Tip	QTY
AEKC1507	7mm	Box of 50
AEKC1508	8mm	Box of 50
AEKC1509	9mm	Box of 50
AEKC1510	10mm	Box of 50
AEKC1511	11mm	Box of 50
AEKC1512	12mm	Box of 50





Pennine Healthcare

Uterine Aspiration Curette System

KARMAN CURETTES

- Clear, flexible curettes for uterine aspiration
- Available with connector to fit to Pennine rotary swivel handle

Product Code	Diameter	QTY
AEAC3004	4mm	Box of 50
AEAC3005	5mm	Box of 50
AEAC3006	6mm	Box of 50
AEAC3007	7mm	Box of 50
AEAC3008	8mm	Box of 50
AEAC3009	9mm	Box of 50
AEAC3010	10mm	Box of 50
AEAC3011	11mm	Box of 50
AEAC3012	12mm	Box of 50

