




## Uterine Aspiration Curette System

### FLEXI-CUFF ASPIRATION TUBING

- Sterile, double wrapped individual peel pouches
- Leak proof seals
- Easily connects to curettes and Karman catheters
- Loose male connector in each pack for connection to suction machines where required
- Made of high density polythene
- Flexi-Cuff provides extra flexibility at a crucial point offering ease of use and improved technique



### Wide Bore Aspiration Tubing (Flexi-Cuff)

Product Number	Description	QTY
KHS-3000	2m wide bore + rotary handle	Each
KHS-3030	3m wide bore + rotary handle	Each



### Wide Bore Aspiration Tubing (without swivel handle and flexi cuff)

Product Number	Description	QTY
AT-2000	2m wide bore	Each
AT-3000	3m wide bore	Each

### ASPIRATION CURETTES


- Clear, rigid plastic curettes for uterine aspiration
- Moulded markings for sizing and correct tip orientation
- Seven sizes available in either curved or straight designs
- Connects to Pennine rotary swivel handle and wide bore, high vacuum tubing

Product Code	Curved Tip	QTY
KC-1007	7mm	Each
KC-1008	8mm	Each
KC-1009	9mm	Each
KC-1010	10mm	Each
KC-1011	11mm	Each
KC-1012	12mm	Each
KC-1014	14mm	Each

Product Code	Straight Tip	QTY
KC-1507	7mm	Each
KC-1508	8mm	Each
KC-1509	9mm	Each
KC-1510	10mm	Each
KC-1511	11mm	Each
KC-1512	12mm	Each
KC-1514	14mm	Each

Product Code	Curved Tip 2m Sets	QTY
KCS-2007	7mm	Box 10
KCS-2008	8mm	Box 10
KCS-2009	9mm	Box 10
KCS-2010	10mm	Box 10
KCS-2011	11mm	Box 10
KCS-2012	12mm	Box 10
KCS-2014	14mm	Box 10

Product Code	Straight Tip 2m Sets	QTY
KCS-2507	7mm	Box 10
KCS-2508	8mm	Box 10
KCS-2509	9mm	Box 10
KCS-2510	10mm	Box 10
KCS-2511	11mm	Box 10
KCS-2512	12mm	Box 10
KCS-2514	14mm	Box 10





## Uterine Aspiration Curette System

### KARMAN CURETTES

- Clear, flexible curettes for uterine aspiration
- Available with connector to fit to Pennine rotary swivel handle

Product Code	Diameter	QTY
AC-3004	4mm	Box of 50
AC-3005	5mm	Box of 50
AC-3006	6mm	Box of 50
AC-3007	7mm	Box of 50
AC-3008	8mm	Box of 50
AC-3009	9mm	Box of 50
AC-3010	10mm	Box of 50
AC-3011	11mm	Box of 50
AC-3012	12mm	Box of 50



### Hy-G-Vac Suction Collection System

A complete range of compatible Pennine products designed for collection of surgical waste during suction procedures

The core of the suction collection system the Pennine Hy-G-Vac collection bottle

- A safe, non-shatter, polypropylene bottle
- Non-removable leak-proof lid
- Quick and easy to assemble and use
- One product for all applications
- Hydrophobic block filter in the lid
- 2 litre capacity
- Safe transport after use including safe incineration properties
- HY-G-Sorb gelling agent to solidify bottle contents before disposal
- Range of trolleys, brackets and other devices to adapt and fix to existing hospital equipment



Product Code	Description	QTY
HGV-2000	Hy-G-Vac suction collection bottle 2 litre capacity	Box 20
HGV-2010	Link (tandem) connecting tube, F/F, 30cm, 6mm	Box 100
HGV-2020	Specimen collection adaptor	Box 10
HGV-2040	Ring holder for bracket or rail attachment	Each
HGV-2050	HY-G-SORB gelling agent (20gm to gel approx. 1000ml water)	Box 20
HGV-2052	HY-G-SORB gelling agent (4.0kg container)	Each