



BladeFLASK - Scalpel Blade Removal System

A simply safer way of removing blades in operating rooms.

Clicksmart BladeFLASK provides a simply safer way to remove and contain used scalpel blades. This world-first single-handed scalpel blade remover prevents the dangerous injuries commonly caused by removing blades with fingers, forceps or re-sheathing.

- Puncture proof sharps container
- Safe to dispose when full
- Holds 100 blades, and has a counter
- Quick, safe, simple.

Remove & Contain Blades with one hand – push and click!

BladeFLASK is secured by a reusable mounting bracket, which means you do not have to hold or brace the BladeFLASK when in use. The bracket enables the BladeFLASK to be positioned at point of use (wall, bench or trolley mounted), making single-handed blade removal, containment and disposal as simple as 'push and click'.

Not just effective... it is cost effective!

This essential safety solution ensures regulatory compliance, dramatically reduces injury costs, and is significantly cheaper than disposable plastic scalpels. You will save up to US\$225 of disposable plastic scalpels with one BladeFLASK.

Sharps container

BladeFLASK is a puncture proof sharps container. When it is full, the BladeFLASK can be disposed safely without any further exposure to sharps contact.



Auto shut-off

BladeFLASK has a built-in counter with a shut-off mechanism that activates automatically after 100 blades have been removed. There is no danger of overfilling.



Product Code	Description	QTY
QFYGEN	Clicksmart BladeFlask Yellow, Blade capacity - 100 blades	Each
OBRGEN	Mounting Bracket for BladeFlask	Each