



Face Masks

Mask With Plastic Dome And Inflatable Cushion - Reusable

LATEX FREE

- With inflatable cushion
- With connector 22mm I.D
- Dome made of hard plastic (PSU)
- Inflatable cushion made of soft silicone
- Offers ideal seal via individual adjustment to the shape of the face



Product Code	Product Description	Size	Colour	QTY
35-70-222	Anaesthetic Face Mask Silicone with Inflatable Air Cushion Cuff Child	2	Green	Each
35-70-223	Anaesthetic Face Mask Silicone with Inflatable Air Cushion Cuff Small Adult	3	Yellow	Each
35-70-224	Anaesthetic Face Mask Silicone with Inflatable Air Cushion Cuff Medium Adult	4	Red	Each
35-70-225	Anaesthetic Face Mask Silicone with Inflatable Air Cushion Cuff Large Adult	5	Purple	Each
35-70-226	Anaesthetic Face Mask Silicone with Inflatable Air Cushion Cuff Extra Large Adult	6	Blue	Each

PVC Face Mask With Inflatable Cushion - Single Use

- With inflatable cushion
- Latex Free



Product Code	Product Description	Size	Colour	QTY
35-40-201	Anaesthesia Face Mask PVC For Newborn, I.D. 15 mm	1	-	Each
35-40-202	Anaesthesia Face Mask PVC For Infant, I.D. 15 mm	2	-	Each
35-40-203	Anaesthesia Face Mask PVC For Child, I.D. 22 mm	3	Yellow	Each
35-40-204	Anaesthesia Face Mask PVC For Small Adult, I.D. 22 mm	4	Red	Each
35-40-205	Anaesthesia Face Mask PVC For Medium Adult, I.D. 22 mm	5	Purple	Each
35-40-206	Anaesthesia Face Mask PVC For Extra Large Adult, I.D. 22 mm	6	Blue	Each

