



# Face Masks

## Bi-Mask With Silicone Lip - Reusable

**LATEX FREE**

- With connector 22mm I.D
- Internal plastic shell (PSU)
- Silicone lip
- Autoclavable



Product Description	Size	Colour	QTY	Product Code
Face Bi-Mask with Silicone Lip Child	2	Green	Each	35-65-222
Face Bi-Mask with Silicone Lip Small Adult	3	Yellow	Each	35-65-223
Face Bi-Mask with Silicone Lip Medium Adult	4	Red	Each	35-65-224
Face Bi-Mask with Silicone Lip Large Adult	5	Purple	Each	35-65-225
Face Bi-Mask with Silicone Lip Extra Large Adult	6	Blue	Each	35-65-226

## PVC Face Mask With Inflatable Cushion - Single Use

- With inflatable cushion
- Latex Free



Product Description	Size	Colour	QTY	Product Code
Face Mask PVC For Newborn, I.D. 15 mm	1	-	Each	35-40-201
Face Mask PVC For Infant, I.D. 15 mm	2	-	Each	35-40-202
Face Mask PVC For Child, I.D. 22 mm	3	Yellow	Each	35-40-203
Face Mask PVC For Small Adult, I.D. 22 mm	4	Red	Each	35-40-204
Face Mask PVC For Medium Adult, I.D. 22 mm	5	Purple	Each	35-40-205
Face Mask PVC For Extra Large Adult, I.D. 22 mm	6	Blue	Each	35-40-206