



# Laryngeal Masks



- The cuff is made of 100% medical grade silicone
- Its special designed shape coincides with the laryngopharynx well
- Reduces stimulation to patient body & improves cuff seal
- Comes in 3 types - Single Use, Reinforced & Reusable
- Latex Free

### Reusable Laryngeal Mask

Autoclave sterilization only. Can be used up to 40 times.

### Reinforced Laryngeal Mask

Spiral reinforcement minimizes crushing or kinking. Smooth transparent and kink resistant tubing.

### Single Use

Product Code	Description	Size	Max Cuff Inflation volume	Qty
LM1004GS1	Neonates up to 5kg	1	4 ml	Each
LM1507GS1	Infants 5-10kg	1.5	7 ml	Each
LM2010GS1	Infants 10-20kg	2	10 ml	Each
LM2514GS1	Children 20-30kg	2.5	14 ml	Each
LM3020GS1	Children 30-50kg	3	20 ml	Each
LM4030GS1	Adults 50-70kg	4	30 ml	Each
LM5040GS1	Adults 70-100kg	5	40 ml	Each
LM1507GS1R	Reinforced Infants 5-10kg	1.5	7 ml	Each
LM2010GS1R	Reinforced Infants 10-20kg	2	10 ml	Each
LM2514GS1R	Reinforced Children 20-30kg	2.5	14 ml	Each
LM3020GS1R	Reinforced Children 30-50kg	3	20 ml	Each
LM4030GS1R	Reinforced Adults 50-70kg	4	30 ml	Each
LM5040GS1R	Reinforced Adults 70-100kg	5	40 ml	Each

### Reusable

Product Code	Description	Size	Max Cuff Inflation volume	Qty
LM1004GS2	Neonates up to 5kg	1	4 ml	Each
LM1507GS2	Infants 5-10kg	1.5	7 ml	Each
LM2010GS2	Infants 10-20kg	2	10 ml	Each
LM2514GS2	Children 20-30kg	2.5	14 ml	Each
LM3020GS2	Children 30-50kg	3	20 ml	Each
LM4030GS2	Adults 50-70kg	4	30 ml	Each
LM5040GS2	Adults 70-100kg	5	40 ml	Each
LM2010GS2R	Reinforced Infants 10-20kg	2	10 ml	Each
LM2514GS2R	Reinforced Children 20-30kg	2.5	14 ml	Each
LM3020GS2R	Reinforced Children 30-50kg	3	20 ml	Each
LM4030GS2R	Reinforced Adults 50-70kg	4	30 ml	Each
LM5040GS2R	Reinforced Adults 70-100kg	5	40 ml	Each

