






Rectoscopes & Insufflator's



Early diagnosis of colorectal cancers and other bowel diseases are made easier with this innovative range of plastic, single-use rectosigmoscopes from Sapi Med. Designed in co-operation with the proctology school of the University of Padova and other leading international institutions, the range includes a choice of different sizes of scopes with tough, hi-tech multi-use plastic grips. Designed for use with the versatile GLF100 cold light sources and fibre-optic cables.

- Range of four sizes available to suit varying clinical needs
- External cm markings on tubes
- Flash-free
- Guaranteed atraumatic
- Cold-light source GLF 100 fits all scopes
- Conform fully with essential requirements of European Council Directive 93/42/EU
- Latex free



Product Code	Description	QTY
A.4113	Short rectoscope, Transparent tube length 130mm, Internal diameter: 17mm, External diameter: 19.6mm	Box 25
A.4125	Rectoscope, Transparent tube, Length: 250mm, Internal diameter: 17mm, External diameter: 20mm	Box 25
A.4130	Rectosigmascope "Rigid-flexible", White tube length: 300mm, Internal diameter: 13.4mm, External diameter: 19.5mm	Box 50
A.4131	Universal Grip	Box 10
A.8200	Disposable Insufflator in PVC Disposable Non Sterile	Each
A.2800ST	Disposable Insufflator in PVC Disposable & Sterile	Each