



# Stylets, Introducer & Tube Exchanger

## I-Bougie

To facilitate orotracheal intubation

- Oxygenate, 2 types of O<sub>2</sub> connectors included
- Centimeter marking from 5 to 40cm
- Unique material:
- Low friction surface
- Allows shaping of the I-Bougie
- For single use, sterile



Size	Length	For ET Tube	QTY	Product Code
14FR	70 cm	≥ 6.0 mm I.D	Each	33-70-700-1

## Tube Exchanger

For extubation and exchange of Endotracheal Tubes

- Lumen-design, O<sub>2</sub> Connector included
- Graduation marks
- For single use, sterile




Size	Length	For ET Tube	QTY	Product Code
11CH	80cm	≥ 4.0 mm	Each	33-45-800-1
14CH	80cm	≥ 5.0 mm	Each	33-46-800-1
19CH	80cm	≥ 7.0 mm	Each	33-47-800-1

## Stylets

Malleable intubation stylet

- Soft, atraumatic tip
- With and without moveable silicone connector
- Reusable



Size	Length	For ET Tube	QTY	Product Code without silicone connector	QTY	Product Code with movable silicone connector (I.D. 15 mm)
Small	34 cm	2.5 - 3.5 mm	Each	90-10-111	Each	90-10-001
Medium	38 cm	4.0 - 4.5 mm	Each	90-10-222	Each	90-10-002
Large	45 cm	5.0 - 6.5 mm	Each	90-10-333	Each	90-10-003
Large	45 cm	7.0 - 11.0 mm	Each	90-10-444	Each	90-10-004

