



Stylets, Introducer & Tube Exchanger

METTS (Muallem ET Tube Stylet)

For difficult intubation

- Core made of metal, malleable (maintains curvature)
- Flexible pre-formed tip
- Graduation marks
- For single use, Sterile

Product Code	Size	Length	For ET Tube	QTY
33-08-400-1	8 CH	40 cm	≥ 3.5 mm	Each
33-12-650-1	12 CH	65 cm	≥ 5.0 mm	Each
33-14-650-1	14 CH	65 cm	≥ 6.0 mm	Each



Pocket Introducer

For difficult intubation

- Folded to 20cm, unfolds to 65cm
- Ready for use within seconds
- Pre-formed tip
- Graduation marks
- For single use, sterile

Product Code	Size	Length	For ET Tube	QTY
33-80-650-1	15 CH	65 cm	≥ 6.0 mm	Each

Pocket Introducer, in robust plastic cartridge

Product Code	Size	Length	For ET Tube	QTY
33-81-650-1	15 CH	65 cm	≥ 6.0 mm	Each



Introducer

For difficult intubation

- Lumen-design, O₂-Connector included
- Rigid pre-formed tip
- Graduation marks
- For single use, Sterile

Product Code	Size	Length	For ET Tube	QTY
33-70-650-1	15 CH	65 cm	≥ 6.0 mm	Each



O2-Connector

For oxygenation, fits all sizes with lumen-design (fits 11 CH, 14 CH, 15 CH, 19 CH), for single use



Product Number	QTY
33-60-005	Box 10





Stylets, Introducer & Tube Exchanger

METTI (Muallem ET Tube Introducer)

- For difficult intubation and tube exchange
- Core made of plastic, semi-rigid (Flexible)
- Flexible, pre-formed tip
- Graduation marks
- For single use, sterile

Product Code	Size	Length	For ET Tube	QTY
33-12-800-1	12 CH	80 cm	≥ 5.0 mm	Each
33-14-800-1	14 CH	80 cm	≥ 6.0 mm	Each



Intubation Catheter

- For fiberoptically assisted intubation (≤ 4.0mm)
- Lumen-design, O₂ Connector included
- Graduation marks
- For single use, sterile

Product Code	Size	Length	For ET Tube	For Fiberoptic	QTY
33-60-560-1	19 CH	56 cm	≥ 7.0 mm	≤ 4.0 mm	Each



Tube Exchanger

- For extubation and exchange of Endotracheal Tubes
- Lumen-design, O₂ Connector included
- Graduation marks
- For single use, sterile

Product Code	Size	Length	For ET Tube	QTY
33-45-800-1	11CH	80cm	≥ 4.0 mm	Each
33-46-800-1	14CH	80cm	≥ 5.0 mm	Each
33-47-800-1	19CH	80cm	≥ 7.0 mm	Each



Stylets

- Malleable intubation stylet
- Soft, atraumatic tip
- With and without moveable silicone connector
- Reusable

Product Code without silicone connector	Size	Length	For ET Tube	QTY	Product Code with movable silicone connector (I.D. 15 mm)	QTY
90-10-111	Small	34 cm	2.5 - 3.5 mm	Each	90-10-001	Each
90-10-222	Medium	38 cm	4.0 - 4.5 mm	Each	90-10-002	Each
90-10-333	Large	45 cm	5.0 - 6.5 mm	Each	90-10-003	Each
90-10-444	Large	45 cm	7.0 - 11.0 mm	Each	90-10-004	Each

