



Suction Coagulator - Precision and Flexible

E Surgical Suction Coagulator is available in both hand switching and foot switching configurations. Suction Coagulator's ergonomic design provides ease of use.

- Hand and foot switching capabilities
- Ergonomically designed for easy handling
- Available in several configurations
- Finger tip suction control

PRODUCT INFORMATION

ACTIVATION	Foot switching or hand switching
SUCTION CONTROL	Regulated by a suction vent on handle
SUCTION TUBE	Aluminium alloy to allow malleable customization
CORD ATTACHED- FOOT SWITCH	One conductor, 10.5 ft length
CORD ATTACHED- HAND SWITCH	Three conductors, 10.5 ft length terminated for direct use with generator



Product Code	Description	QTY
SC08F	Foot Switching 13 cm (5 in) 8 FR	Box 25
SC10F	Foot Switching 13 cm (5 in) 10 FR	Box 25
SC12F	Foot Switching 13 cm (5 in) 12 FR	Box 25
SC08H	Hand Switching 13 cm (5 in) 8 FR	Box 25
SC10H	Hand Switching 13 cm (5 in) 10 FR	Box 25
SC12H	Hand Switching 13 cm (5 in) 12 FR	Box 25
AT300	Tip Cleaner 42.5mm x 42.5mm	Box 50
AH300	Holster Standard Size	Box 50