




Endotracheal Tubes

TRACHEAL TUBE WITHOUT CUFF

- Transparent thermosensitive PVC
- Oral/nasal
- Murphy eye
- X-ray line
- Supplied sterile, single-use, boxes of 10



ID (mm)	OD (mm)	Code PVC
2.0	3.0	22-2000
2.5	3.5	22-2500
3.0	4.2	22-3000
3.5	4.8	22-3500
4.0	5.4	22-4000
4.5	6.0	22-4500
5.0	6.9	22-5000
5.5	7.5	22-5500
6.0	8.2	22-6000

ID (mm)	OD (mm)	Code PVC
6.5	8.8	22-6500
7.0	9.7	22-7000
7.5	10.3	22-7500
8.0	11.0	22-8000
8.5	11.7	22-8500
9.0	12.3	22-9000
9.5	13.0	22-9500
10.0	13.7	22-1000

TRACHEAL TUBE WITH HIGH VOLUME CUFF

- Transparent, thermosensitive PVC
- Oral/nasal
- Thin-walled high volume cuff and marked pilot balloon
- Ideal for long-term ventilation
- Murphy eye
- X-ray line
- Supplied sterile, single-use, boxes of 10



ID (mm)	OD (mm)	Cuff (mm)	Code
3.0	4.2	22	24-3095
3.5	4.9	22	24-3595
4.0	5.5	22	24-4095
4.5	6.2	22	24-4595
5.0	6.9	25	24-5095
5.5	7.5	25	24-5595
6.0	8.2	30	24-6095
6.5	8.8	30	24-6595
7.0	9.7	34	24-7095
7.5	10.3	34	24-7595

ID (mm)	OD (mm)	Cuff (mm)	Code
8.0	11.0	34	24-8095
8.5	11.7	34	24-8595
9.0	12.3	35	24-9095
9.5	13.0	35	24-9595
10.0	13.7	35	24-1095



Endotracheal Tubes

REINFORCED TRACHEAL TUBE

- Oral/nasal
- With stainless steel kink-resistant reinforcement
- Murphy eye
- With marked pilot balloon
- Supplied sterile, single-use, boxes of 10



ID (mm)	OD (mm)	Cuff (mm)	Code with cuff PVC	Code without cuff PVC
3.0	5.3	8	33-3000	26-3000
3.5	5.8	10	33-3500	26-3500
4.0	6.3	10	33-4000	26-4000
4.5	6.8	12	33-4500	26-4500
5.0	7.3	14	33-5000	26-5000
5.5	7.8	17	33-5500	26-5500
6.0	8.4	20	33-6000	26-6000
6.5	9.6	20	33-6500	26-6500

ID (mm)	OD (mm)	Cuff (mm)	Code with Cuff PVC	Code without Cuff PVC
7.0	10.2	23	33-7000	26-7000
7.5	11.0	25	33-7500	26-7500
8.0	11.6	26	33-8000	26-8000
8.5	12.2	27	33-8500	26-8500
9.0	12.9	29	33-9000	26-9000

TRACHEAL TUBE WITH SUCTION LUMEN

- Allows suction of secretion from above the cuff, minimising risk of infection
- Silicone PVC
- Oral/nasal
- Delicate low-pressure cuff
- Murphy eye
- Available with cap or connector
- With marked pilot balloon and x-ray line
- Supplied sterile, single-use, boxes of 10



ID (mm)	OD (mm)	Cuff (mm)	Code with cap
6.0	8.2	20	11-6015
6.5	8.8	20	11-6515
7.0	9.7	23	11-7015
7.5	10.3	25	11-7515


ID (mm)	OD (mm)	Cuff (mm)	Code with cap
8.0	11.0	26	11-8015
8.5	11.7	27	11-8515
9.0	12.3	29	11-9015



Endotracheal Tubes

PREFORMED TRACHEAL TUBE, NASAL, NORTH

- Soft flexible PVC
- Murphy eye
- X-ray line and marked pilot balloon
- Supplied sterile, single-use, boxes of 10



ID (mm)	OD (mm)	Cuff (mm)	Code with cuff PVC	Code without cuff PVC
3.0	4.2	8	38-3000	40-3000
3.5	4.8	10	38-3500	40-3500
4.0	5.4	10	38-4000	40-4000
4.5	6.0	12	38-4500	40-4500
5.0	6.9	14	38-5000	40-5000
5.5	7.5	17	38-5500	40-5500

ID (mm)	OD (mm)	Cuff (mm)	Code with cuff PVC	Code without cuff PVC
6.0	8.2	20	38-6000	40-6000
6.5	8.8	20	38-6500	40-6500
7.0	9.7	23	38-7000	40-7000
7.5	10.3	25	38-7500	40-7500
8.0	11.0	26	38-8000	40-8000

PREFORMED TRACHEAL TUBE, ORAL, SOUTH

- Available with or without low pressure cuff
- Murphy eye
- With x-ray line and marked pilot balloon
- Supplied sterile, single-use, boxes of 10



ID (mm)	OD (mm)	Cuff (mm)	Code with cuff PVC	Code without cuff PVC
3.0	4.2	8	34-3000	37-3000
3.5	4.8	10	34-3500	37-3500
4.0	5.4	10	34-4000	37-4000
4.5	6.0	12	34-4500	37-4500
5.0	6.9	14	34-5000	37-5000
5.5	7.5	17	34-5500	37-5500
6.0	8.2	20	34-6000	37-6000

ID (mm)	OD (mm)	Cuff (mm)	Code with cuff PVC	Code without cuff PVC
6.5	8.8	20	34-6500	37-6500
7.0	9.7	23	34-7000	37-7000
7.5	10.3	25	34-7500	37-7500
8.0	11.0	26	34-8000	37-8000
8.5	11.7	27	34-8500	37-8500
9.0	12.3	29	34-9000	37-9000