



Oxygen Delivery Kink Resistant Tubing

Oxygen tubes are used to connect the oxygen source with the oxygenation device during administration of oxygen. The product is intended for single use only. Innovative 3-channel safety tubing is designed to maintain an uninterrupted flow of oxygen, including our patented tubing end fitting.



2025



- Salter's exclusive 3-channel safety tubing reduces the risk of oxygen flow loss due to twisting and kinking
- Soft, stretchable tubing lessens pressure on the ears and face
- Feels light to the touch
- Non-DEHP throughout
- Redesigned slide bolo helps resist twisting
- Maximizes comfort and wear-ability

Line Length	QTY	Product Code Clear	Product Code Green
7' (2.1m)	Box 50	2002	2002G
9' (2.7m)	Box 50	2009	
10' (3m)	Box 50	2010	2010G
15' (4.5m)	Box 50		2015G
20' (6m)	Box 50	2020	2020G
25' (7.3m)	Box 50	2025	2025G
30' (9.1m)	Box 50	2030	
35' (10.6m)	Box 50	2035	2035G
40' (12.1m)	Box 50	2040	2040G
50' (15.2m)	Box 50	2050	2050G

