



Pelvic Sling

PNEUMATIC SYSTEM FOR STABILISATION OF PELVIC FRACTURES

The pneumatic Pelvic Sling provides early stabilisation of pelvic fractures at the accident scene to allow patient transport and to reduce life threatening bleeding.

The inflatable cuffs exert a lateral pressure to the pelvis for external stabilisation of the pelvic ring with the purpose to regain the anatomic conditions. This reduces the inner volume and tamponades the bleeding. The integrated inflatable cuff exerts a pressure from 70 to 100 mmHg and helps to reduce bleeding and pain significantly.

Slide Spint / Knee Fixation

to slide the Pelvic Sling underneath the patients pelvis. After removal used as knee fixation.



Cuffs

Both lateral cuffs are connected and can be inflated with hand inflator for an ideal stabilisation of the pelvis. No force required due to inflatable cuffs.

Product Code	Description	QTY
34-10-006	Pelvic Sling with manometer and hand inflator, size S junior, pelvic circumference 70-90cm, single use	Each
34-10-004	Pelvic Sling with manometer and hand inflator, size M standard, pelvic circumference 90-110cm, single use	Each
34-10-005	Pelvic Sling with manometer and hand inflator, size L large, pelvic circumference 110-140cm, single use	Each