



Airway Management Simulators

By nature airway management includes uncertainty, complexity and time pressure. Simulation is the ideal method for replicating clinical situation and contributes to improve patient care by enhancing proficiency of practitioners. VBM has developed a wide range of skill-trainers and mannequins with the highest possible fidelity to train all airway procedures. The devices cover clinical situations such as direct laryngoscopy, video laryngoscopy, fiberoptic bronchoscopy, supraglottic airway devices insertion, Trans Tracheal Jet-Ventilation, bagmask ventilation and cricothyrotomy.

Especially cricothyrotomy is a rarely performed procedure the technical proficiency can realistically only be obtained through simulated learning. VBM offers the widest range of cricothyrotomy skill-trainers for simulation.



Product Number	Description	QTY
30-29-000	BILL - ET Intubation - Supraglottic airway insertion - Mask ventilation - Fiberoptic bronchoscopy (if bronchial tree is connected) - Fiberoptic intubation - Video laryngoscopy - Inflatable tongue for Difficult Airway simulation Complete with base, carina with bayonet lock and carrying bag	Each



Product Number	Description	QTY
30-19-400	Bronchial tree acc. to Dr. Nakhosteen For connection to simulator „BILL“	Each



Product Number	Description	QTY
30-16-300	Transparent Head Cross section with upper airways For demonstration and positioning of Laryngeal Tube # 3, 4 and iLTS-D # 4/5, in carrying bag	Each



Airway Management Simulators



Product Code	Description	QTY
30-14-700	Crico-Trainer "Adelaide" The new system allows training of surgical cricothyrotomy techniques such as open and wire-guided cricothyrotomy. The special concept of skin with 2 layers is simulating the subcutaneous tissue and the cricothyroid membrane. It is designed to replicate difficult scenario thanks to the moveable chin. Complete with 5 skins	Each
30-14-000	Crico-Trainer "Frova" For cricothyrotomy and percutaneous tracheostomy, Complete with 10 skins and 10 tracheas	Each
30-14-900	Crico-Trainer "Paediatric" Complete with 10 skins	Each
30-14-500	Crico-Trainer "Pig" for fixation of a real animal trachea Complete with 10 skins	Each
30-14-222	Skin For Crico-Trainer "Frova" and "Pig"	Each
30-14-722	Skin For Crico-Trainer "Adelaide"	Each
30-14-922	Skin For Crico-Trainer "Paediatric"	Box 10
30-14-111	Larynx For Crico-Trainer "Frova" and "Adelaide"	Each
30-14-444	Trachea For Crico-Trainer "Frova"	Each