



# Link Yankauers

## Link Yankauers

- Sterile, double wrapped transparent individual peel pouches
- Available with or without vacuum control (VC) vent
- Short mini, Mini fine, Midi and Maxi
- Yankauers available with smooth rounded end or Super Tip
- Midi Yankauer available with additional Atrauma Tip to reduce tissue grab

Product Code	Type	Length	CH	Use	QTY
LY-3400	Short Mini	190mm	15CH	Paediatric/Dental	Each
LY-3401	Short Mini VC	190mm	15CH	Paediatric/Dental	Each
LY-4400	Short Mini Super-tip	195mm	15CH	Paediatric/Dental	Each
LY-4401	Short Mini Super-tip VC	195mm	15CH	Paediatric/Dental	Each
LY-3500	Mini Fine Tip	270mm	12CH	Fine Suctioning (e.g ENT)	Each
LY-3501	Mini Fine Tip VC	270mm	12CH	Fine Suctioning (e.g ENT)	Each
LY-4500	Mini Fine Super-tip	275mm	12CH	Fine Suctioning (e.g ENT)	Each
LY-4501	Mini Fine Super-tip VC	275mm	12CH	Fine Suctioning (e.g ENT)	Each
LY-3600	Midi 30o angle	280mm	18CH	For general and pharyngeal suction	Each
LY-3601	Midi VC	280mm	18CH	For general and pharyngeal suction	Each
LY-3600/PR	Midi 10o angle	270mm	18CH	For general and pharyngeal suction	Each
LY-4600	Midi Super-tip	280mm	18CH	For general and pharyngeal suction	Each
LY-4601	Midi Super-tip VC	280mm	18CH	For general and pharyngeal suction	Each
LY-3700	Maxi	275mm	22CH	Large particles and significant volumes of fluid	Each
LY-3701	Maxi VC	275mm	22CH	Large particles and significant volumes of fluid	Each
LY-4700	Maxi Super-tip	275mm	22CH	Large particles and significant volumes of fluid	Each
LY-4701	Maxi Super-tip VC	275mm	22CH	Large particles and significant volumes of fluid	Each
LY-3800	Maxi Plus	280mm	28CH	For orthopaedic use	Each
LY-3801	Maxi Plus VC	280mm	28CH	For orthopaedic use	Each
AELY3602	Atrauma tip, short soft	280mm	18CH	General suction with minimal tissue grab	Each
AELY3603	As above with VC	280mm	18CH	General suction with minimal tissue grab	Each
AELY3606	Atrauma tip, soft rigid	280mm	18CH	General suction with minimal tissue grab	Each
AELY3607	As above with VC	280mm	18CH	General suction with minimal tissue grab	Each
AELY3608	Atrauma tip, long rigid	280mm	18CH	General suction with minimal tissue grab	Each
AELY3609	As above with VC	280mm	18CH	General suction with minimal tissue grab	Each

