



Bile Bags



LATEX FREE

~~PHT
DEHP~~
DEHP-free

- For collection of bile and other gastric secretions
- Easy connection to a Pennine nasogastric tube (Ryles / Levins tube) or with Kehr's tube after bile duct surgery
- Stop cock for cleaner emptying procedure
- Available with latex straps for bed or patient attachment
- All angles, lengths & dimensions are nominal

Product Code	Description	QTY
AEBB350RY	Ryles / Levin Tube connections with stop cock	Each
AEBB350TR	Naso Gastric Drainage Bag, Ryles Connector, & T Connector	Each
AEBB350TY	Kehr's tube or Female Luer Connection with stop cock	Each

