



Feeding Tubes

Graduated Infant Feeding Tubes

- Tube with radio-opaque line provided to its full length
- Sterile, individual peel packs
- Single use
- Closed tip, 2 side eyes to reduce risk of blockage
- Color coded funnel and connector for easy identification of size



Graduated X-ray tube

CH/Dia.	QTY	Product Code 40cm	Product Code 50m	Product Code 100cm
5ch ø 1.7mm	Each	AEFT190540	AEFT190550	AEFT1905100
6ch ø 2.00mm	Each	AEFT190640	AEFT190650	AEFT1906100
8ch ø 2.7mm	Each	AEFT190840	AEFT190850	AEFT1908100
10ch ø 3.3mm	Each	AEFT191040	AEFT191050	