



Floor Suction Devices

The Quiet Floor Suction range is a collection of floor suction solutions that have a dual action – double effect by using both superior wicking and vacuum suction to achieve low noise and a drier operating room floor.



WaterBug Quiet Floor Suction® Device

- Efficiently suctions fluid waste to help reduce risks and maintain a dry operating room floor
- Its raised design and 365cm of flexible tubing allows for easy repositioning
- Suctions unlimited amounts of irrigants, blood, and other bodily fluids from the floor
- For thorough pickup, fluid is captured from the top via holes as well as from the bottom through suction and wicking action
- Non-Sterile, intended for single-use, and not made with natural

Product Code	Description	QTY
90030	WaterBug Quiet Floor Suction® Device With Tube	Each



WaterBoom® Floor Suction Strip

- Creates a dam to help hold back, contain, and quietly suction an unlimited amount of waste fluid
- Provides strong suction, while conforming to your O.R. floor
- Its low profile helps you step over it without suction loss
- Its 365cm tubing attaches from the side to help reduce tripping hazards
- This unobtrusive device does not require repositioning and is a good choice for those who prefer nothing at their feet
- Non-Sterile, intended for single-use, and not made with natural rubber latex

Product Code	Description	QTY
94030	WaterBoom® Floor Suction Strip 81cm With Tube	Each