




Laparoscopy For Gynaecology

Suction Irrigation Set

- Transparent ergonomic handpiece - comfort design for both right and left hand use
- Colour coded smooth push button pistons with dual seal to prevent leaking
- High flow clog free valves
- 5mmx 33cm rigid probe
- Dual spike with 3 metres of flexible tubing
- Flip back cap (fluid and gas tight) facilitates the utilisation of instruments up to 3.5mm or laser fibres through the hand-piece.
- Smooth atraumatic probe tip with 8 vents
- Quick probe exchange (10mm probe available LAP0003)
- Latex Free



Product Code	Description	QTY
LAP0008A	Suction irrigation set - 5mm probe	Box 10
SI0006	10mm probe for suction irrigation set	Box 10