




# Levins/Duodenal Tubes

## Levins/Duodenal Tubes

Duodenal tubes are for active aspiration, feeding, injection, lavage or other gastric surgical procedures. With or without X-Ray detectable line.

- Soft, frosted and kink resistance PVC tubing
- Sterile, individual peel pouches
- 4 side eyes to reduce risk of blockage
- Open terminal eye to ensure maximum dispersion
- Available with X-Ray opaque line
- Latex free



Levins Tubes - X-Ray line

| Open end & X-Ray line | QTY  | CH/Diam      |
|-----------------------|------|--------------|
| AEXLT2206             | Each | 6ch ø 2.0mm  |
| AEXLT2208             | Each | 8ch ø 2.7mm  |
| AEXLT2210             | Each | 10ch ø 3.3mm |
| AEXLT2212             | Each | 12ch ø 4.0mm |
| AEXLT2214             | Each | 14ch ø 4.7mm |
| AEXLT2216             | Each | 16ch ø 5.3mm |
| AEXLT2218             | Each | 18ch ø 6.0mm |
| AEXLT2220             | Each | 20ch ø 6.7mm |