




Levins/Duodenal Tubes


Levins/Duodenal Tubes

Duodenal tubes are for active aspiration, feeding, injection, lavage or other gastric surgical procedures. With or without X-Ray detectable line.

- Medical grade PVC (Shore Hardness A80)
- Sterile, individual peel pouches
- Printed indications at 40, 45, 50, 55, 60, 65 & 70 cm from tip
- 4 side eyes to reduce risk of blockage
- Terminal eye to ensure maximum dispersion
- Optional closed end
- Available with X-Ray opaque line, XLT type
- Latex free



Closed end	Open end & X-Ray line	Closed end & X-Ray line	QTY	CH/Diam
LT-2406	XLT-2206	XLT-2406	Each	6ch ø 2.00mm
LT-2408	XLT-2208	XLT-2408	Each	8ch ø 2.66mm
LT-2410	XLT-2210	XLT-2410	Each	10ch ø 3.33mm
LT-2412	XLT-2212	XLT-2412	Each	12ch ø 4.00mm
LT-2414	XLT-2214	XLT-2414	Each	14ch ø 4.66mm
LT-2416	XLT-2216	XLT-2416	Each	16ch ø 5.33mm
LT-2418	XLT-2218	XLT-2418	Each	18ch ø 6.00mm
LT-2420	XLT-2220	XLT-2420	Each	20ch ø 6.66mm



Levins Feeding Nasogastric Tube

X-Ray Graduated with ENfit connector, open tip and four side eyes

- Printed indications at 40-45-50-55-60-65 and 70cm from tip
- Terminal eye to ensure maximum dispersion

Open end & X-Ray line	QTY	CH/Diam
XLT-2205/105/EF	Each	5ch ø 1.66mm
XLT-2206/EF	Each	6ch ø 2.00mm
XLT-2208/EF	Each	8ch ø 2.66mm
XLT-2210/EF	Each	10ch ø 3.33mm
XLT-2212/EF	Each	12ch ø 4.00mm
XLT-2214/EF	Each	14ch ø 4.66mm
XLT-2216/EF	Each	16ch ø 5.33mm
XLT-2218/EF	Each	18ch ø 6.00mm
XLT-2220/EF	Each	20ch ø 6.66mm

New product to Bosco - 8 weeks to order in