



Link Yankauers

Link Yankauers

- Sterile, double wrapped transparent individual peel pouches
- Available with or without vacuum control (VC) vent
- Short mini, Mini fine, Midi and Maxi
- Yankauers available with smooth rounded end or Super Tip
- Midi Yankauer available with additional Atrauma Tip to reduce tissue grab

Product Code	Type	Length	CH	Use	QTY
LY-3400	Short Mini	190mm	15CH	Paediatric/Dental	Each
LY-3401	Short Mini VC	190mm	15CH	Paediatric/Dental	Each
LY-4400	Short Mini Super-tip	195mm	15CH	Paediatric/Dental	Each
LY-4401	Short Mini Super-tip VC	195mm	15CH	Paediatric/Dental	Each
LY-3500	Mini Fine Tip	270mm	12CH	Fine Suctioning (e.g ENT)	Each
LY-3501	Mini Fine Tip VC	270mm	12CH	Fine Suctioning (e.g ENT)	Each
LY-4500	Mini Fine Super-tip	275mm	12CH	Fine Suctioning (e.g ENT)	Each
LY-4501	Mini Fine Super-tip VC	275mm	12CH	Fine Suctioning (e.g ENT)	Each
LY-3600	Midi 30o angle	280mm	18CH	For general and pharyngeal suction	Each
LY-3601	Midi VC	280mm	18CH	For general and pharyngeal suction	Each
LY-3600/PR	Midi 10o angle	270mm	18CH	For general and pharyngeal suction	Each
LY-4600	Midi Super-tip	280mm	18CH	For general and pharyngeal suction	Each
LY-4601	Midi Super-tip VC	280mm	18CH	For general and pharyngeal suction	Each
LY-3700	Maxi	275mm	22CH	Large particles and significant volumes of fluid	Each
LY-3701	Maxi VC	275mm	22CH	Large particles and significant volumes of fluid	Each
LY-4700	Maxi Super-tip	275mm	22CH	Large particles and significant volumes of fluid	Each
LY-4701	Maxi Super-tip VC	275mm	22CH	Large particles and significant volumes of fluid	Each
LY-3800	Maxi Plus	280mm	28CH	For orthopaedic use	Each
LY-3801	Maxi Plus VC	280mm	28CH	For orthopaedic use	Each
LY-3602	Atrauma tip, short soft	280mm	18CH	General suction with minimal tissue grab	Each
LY-3603	As above with VC	280mm	18CH	General suction with minimal tissue grab	Each
LY-3606	Atrauma tip, soft rigid	280mm	18CH	General suction with minimal tissue grab	Each
LY-3607	As above with VC	280mm	18CH	General suction with minimal tissue grab	Each
LY-3608	Atrauma tip, long rigid	280mm	18CH	General suction with minimal tissue grab	Each
LY-3609	As above with VC	280mm	18CH	General suction with minimal tissue grab	Each

