



# MPR (Multi-Purpose Rectoscopes)

## Multi-Purpose Rectoscopes for endorectal ecography and biopsies

The main characteristic of MPR rectoscopes is their grip. It features a special opening which allows the insertion of biopsy forceps avoiding the air insufflated to come out. Set on the grip window is a special transparent plastic component allowing featuring a 2x zoom.

There are two versions:

MPR: 220 mm long;  
18 mm of internal diameter

MPR Short: 190 mm long;  
18 mm of internal diameter



| Product Code | Description                        | QTY  |
|--------------|------------------------------------|------|
| A.4519       | MPR Short Multi-purpose Rectoscope | Each |
| A.4522       | MPR Long Multi-purpose Rectoscope  | Each |