



Oppo-Cath Suction Catheters

Oppo-Cath Suction Catheters

Suction catheters for oral, tracheal and endobronchial suction procedures. For increased patient safety and reduced tissue damage. To be used by connecting to the suction source with Pennine suction connecting tubing, FF type



- 48cm/60cm
- PVC Medical grade
- Clear Prestrol Vacuum control
- Graduations every 1cm
- Sterile, individual peel pouch
- Terminal eye, two small opposed side eyes

Product Code 48cm	Product Code 60cm	CH/Diam	QTY
OSC-1205/48	OSC-1205/60	5ch ø 1.7mm	Each
OSC-1206/48	OSC-1206/60	6ch ø 2.0mm	Each
OSC-1207/48	OSC-1207/60	7ch ø 2.3mm	Each
OSC-1208/48	OSC-1208/60	8ch ø 2.7mm	Each
OSC-1210/48	OSC-1210/60	10ch ø 3.3mm	Each
OSC-1212/48	OSC-1212/60	12ch ø 4.00mm	Each
OSC-1214/48	OSC-1214/60	14ch ø 4.7mm	Each
OSC-1216/48	OSC-1216/60	16ch ø 5.3mm	Each

Prestrol Oppo-Cath Suction Catheters

Suction catheters for oral, tracheal and endobronchial suction procedures. For increased patient safety and reduced tissue damage. To be used by connecting to the suction source with Pennine suction connecting tubing, FF type

- PVC Medical grade
- Terminal eye, two small opposed side eyes
- Sterile, individual peel pouches
- With colour coded Prestrol Vacuum Control device

Product Code 48cm	Product Code 60cm	CH/Diam	QTY
OSC-1505/48	OSC-1505/60	5ch ø 1.7mm	Each
OSC-1506/48	OSC-1506/60	6ch ø 2.0mm	Each
OSC-1508/48	OSC-1508/60	8ch ø 2.7mm	Each
OSC-1510/48	OSC-1510/60	10ch ø 3.3mm	Each
OSC-1512/48	OSC-1512/60	12ch ø 4.0mm	Each
OSC-1514/48	OSC-1514/60	14ch ø 4.7mm	Each
OSC-1516/48	OSC-1516/60	16ch ø 5.3mm	Each
OSC-1518/48	OSC-1518/60	18ch ø 6.0mm	Each