



# Oppo-Cath Suction Catheters

## Oppo-Cath Suction Catheters

Suction catheters for oral, tracheal and endobronchial suction procedures. For increased patient safety and reduced tissue damage. To be used by connecting to the suction source with Pennine suction connecting tubing, FF type

- 48cm/60cm
- PVC Medical grade
- Clear Prestrol Vacuum control
- Graduations every 1cm
- Sterile, individual peel pouch
- Terminal eye, two small opposed side eyes



Product Code 48cm	Product Code 60cm	CH/Diam	QTY
OSC-1205/48	-	5ch ø 1.66mm	Each
OSC-1206/48	-	10ch ø 3.33mm	Each
OSC-1207/48	-	12ch ø 4.00mm	Each
OSC-1208/48	OSC-1208/60	8ch ø 2.66mm	Each
OSC-1210/48	OSC-1210/60	10ch ø 3.33mm	Each
OSC-1212/48	OSC-1212/60	12ch ø 4.00mm	Each
OSC-1214/48	OSC-1214/60	14ch ø 4.66mm	Each
OSC-1216/48	OSC-1216/60	16ch ø 5.33mm	Each

## Prestrol Oppo-Cath Suction Catheters

Suction catheters for oral, tracheal and endobronchial suction procedures. For increased patient safety and reduced tissue damage. To be used by connecting to the suction source with Pennine suction connecting tubing, FF type

- PVC Medical grade
- Terminal eye, two small opposed side eyes
- Sterile, individual peel pouches
- With colour coded Prestrol Vacuum Control device



## Standard PVC

Product Code 48cm	Product Code 60cm	CH/Diam	QTY
OSC-1505/48	-	5ch ø 1.66mm	Each
OSC-1506/48	-	6ch ø 2.00mm	Each
OSC-1508/48	OSC-1508/60	8ch ø 2.66mm	Each

## Soft PVC

Product Code 48cm	Product Code 60cm	CH/Diam	QTY
OSC-1510/48	OSC-1510/60	10ch ø 3.33mm	Each
OSC-1512/48	OSC-1512/60	12ch ø 4.00mm	Each
OSC-1514/48	OSC-1514/60	14ch ø 4.66mm	Each
OSC-1516/48	OSC-1516/60	16ch ø 5.33mm	Each
OSC-1518/48	OSC-1518/60	18ch ø 6.00mm	Each