



Pool Suckers & Selectovac Yankauer

Pool Suckers

Flexible or Rigid sump sucker with removable inner suction probe. Specifically for abdominal surgery, reducing trauma to surrounding tissues.

- Available with atraumatic end
- Available with angled tip
- Removable inner probe

Product Code	Type	Length	CH/Tip	Inside Diam.	QTY
YS-3006	Flexible straight tip	300mm	20CH	6.66mm	Each
AEYS3007	Rigid (White)	300mm	15CH	6.5mm	Each
YS-3008	Flexible (angled tip)	300mm	20CH	6.66mm	Each



Flexible straight tip



Flexible angled tip



Rigid tip

Selectovac Yankauer

Surgical suction cannula with switch on and off valve, available with bulb tip and straight tip.

- Switch system allows the operator to manually open and close suction
- Easy connection to Pennine range of tubing



Product Code	Type	Length	CH/Tip	Inside Diam.	QTY
AEYS3009	Bulb tip	265mm	24CH	8mm	Each
AEYS3010	Straight Tip	265mm	24CH	8mm	Each