



## Pressure Infusion Cuff - Reusable

### INFUSOR (500 ml / 1000 ml / 3000 ml)

- Allows clear view on the fluid bag and the fluid level
- 360° visible
- Fluid bag can be replaced quickly and easily
- Consistent pressure on the complete fluid bag
- Precision manometer
- Precise pressure range from 0-300 mmHg
- Large volume allows quick inflation
- For pressure control
- For quick deflation
- Easy to clean (PU material)
- Latex Free



Size	QTY	Product Code With Female Luer Lock	Product Code With Hand Inflator, Stop Cock And Colour Coded Pressure Indicator	Product Code With Hand Inflator And Precision Manometer
500ml	Each	-	57-02-050	57-08-050
1000ml	Each	-	57-02-100	57-08-100
3000ml	Each	57-01-300		57-08-300

### PRESSURE INFUSION CUFFS 1000 ml

- Allows clear view on the fluid bag and the fluid level
- Fluid bag can be replaced quickly and easily
- Consistent pressure on the complete fluid bag
- Precision manometer
- Precise pressure range from 0-300 mmHg
- Large volume allows quick inflation
- With release valve for pressure control



Size	QTY	Product Code With Hand Inflator, Stop Cock And Colour Coded Pressure Indicator	Product Code With Hand Inflator And Precision Manometer
1000ml	Each	58-02-100	58-08-100

### PRESSURE INFUSION CUFFS 500 ml

- Allows clear view on the fluid bag and the fluid level
- Fluid bag can be replaced quickly and easily
- Precision manometer
- Precise pressure range from 0-300 mmHg
- Large volume allows quick inflation
- With release valve for pressure control



Size	QTY	Product Code With Hand Inflator, Stop Cock And Colour Coded Pressure Indicator	Product Code With Hand Inflator And Precision Manometer
500ml	Each	58-02-050	58-08-050



## Pressure Infusion Cuff - Reusable

### PRESSURE INFUSION CUFFS 500 ml

- Allows clear view on the fluid bag and the fluid level
- Fluid bag can be replaced quickly and easily
- Precision manometer
- Precise pressure range from 0-300 mmHg
- Large volume allows quick inflation
- With release valve for pressure control



Size	QTY	Product Code With Hand Inflator, Stop Cock And Colour Cod-ed Pressure Indicator	Product Code With Hand Inflator And Precision Manometer
500ml	Each	58-02-050	58-08-050

### PRESSURE INFUSION CUFFS 3000 / 5000 ml

- Allows clear view on the fluid bag and the fluid level
- Pressure Infusion Cuffs for irrigation purposes during minimal invasive surgery
- Pocket design with zip
- Precise pressure range from 0-300 mmHg
- Large volume allows quick inflation
- With release valve for pressure control



Size	QTY	Product Code With Hand Inflator, Stop Cock And Colour Cod-ed Pressure Indicator	Product Code With Hand Inflator And Precision Manometer
1000ml	Each	58-02-100	58-08-100



## Pressure Infusion Cuff - Disposable

**Disposable (500 ml / 1000 ml / 3000 ml)**

**Transparent Netting**

- Allows clear view on the fluid bag and the fluid level

**Pocket Design with hanger/hook\***

- Fluid bag can be replaced quickly and easily

**Pressure Gauge**

- Colour coded with overpressure valve
- 360° visible

**Hand Inflator**

- Large volume allows quick inflation

**Three-Way Stop Cock**

- For pressure control
- For quick deflation

**For single use**

- Reduces the potential for cross-contamination
- Latex Free



Size	QTY	Product Code With Female Luer Lock	Product Code With Hand Inflator, Stop Cock And Colour Coded Pressure Indicator
500ml	Box 5		56-02-050
1000ml	Box 5	56-01-100	56-02-100
3000ml	Each	56-01-300	56-02-300