



## Pressure Infusion Cuff - Reusable

### INFUSOR (500 ml / 1000 ml / 3000 ml)

- Allows clear view on the fluid bag and the fluid level
- 360° visible
- Fluid bag can be replaced quickly and easily
- Consistent pressure on the complete fluid bag
- Precision manometer
- Precise pressure range from 0-300 mmHg
- Large volume allows quick inflation
- For pressure control
- For quick deflation
- Easy to clean (PU material)
- Latex Free



Product Code With Female Luer Lock	Product Code With Hand Inflator, Stop Cock And Colour Coded Pressure Indicator	Product Code With Hand Inflator And Precision Manometer	Size	QTY
-	57-02-050	57-08-050	500ml	Each
-	57-02-100	57-08-100	1000ml	Each
57-01-300	57-02-300	57-08-300	3000ml	Each

### PRESSURE INFUSION CUFFS 1000 ml

- Allows clear view on the fluid bag and the fluid level
- Fluid bag can be replaced quickly and easily
- Consistent pressure on the complete fluid bag
- Precision manometer
- Precise pressure range from 0-300 mmHg
- Large volume allows quick inflation
- With release valve for pressure control



Product Code With Hand Inflator, Stop Cock And Colour Coded Pressure Indicator	Product Code With Hand Inflator And Precision Manometer	Size	QTY
58-02-100	58-08-100	1000ml	Each

### PRESSURE INFUSION CUFFS 500 ml

- Allows clear view on the fluid bag and the fluid level
- Fluid bag can be replaced quickly and easily
- Precision manometer
- Precise pressure range from 0-300 mmHg
- Large volume allows quick inflation
- With release valve for pressure control



Product Code With Hand Inflator, Stop Cock And Colour Coded Pressure Indicator	Product Code With Hand Inflator And Precision Manometer	Size	QTY
58-02-050	58-08-050	500ml	Each



## Pressure Infusion Cuff - Reusable

### PRESSURE INFUSION CUFFS 500 ml

- Allows clear view on the fluid bag and the fluid level
- Fluid bag can be replaced quickly and easily
- Precision manometer
- Precise pressure range from 0-300 mmHg
- Large volume allows quick inflation
- With release valve for pressure control

Product Code With Hand Inflator, Stop Cock And Colour Coded Pressure Indicator	Product Code With Hand Inflator And Precision Manometer	Size	QTY
58-02-050	58-08-050	500ml	Each



### PRESSURE INFUSION CUFFS 1500 ml

- Allows clear view on the fluid bag and the fluid level
- Pocket design with hook
- Precision manometer
- Precise pressure range from 0-300 mmHg
- Large volume allows quick inflation
- With release valve for pressure control
- Bladder is replaceable

With Female Luer Lock	Product Code With Hand Inflator And Precision Manometer	Size	QTY
53-01-150	53-08-150	1500ml	Each



### PRESSURE INFUSION CUFFS 3000 / 5000 ml

- Allows clear view on the fluid bag and the fluid level
- Pressure Infusion Cuffs for irrigation purposes during minimal invasive surgery
- Pocket design with zip
- Precise pressure range from 0-300 mmHg
- Large volume allows quick inflation
- With release valve for pressure control
- Bladder is replaceable

Product Code With Hand Inflator, Stop Cock And Colour Coded Pressure Indicator	Product Code With Hand Inflator And Precision Manometer	Size	QTY
58-02-100	58-08-100	1000ml	Each





## Pressure Infusion Cuff - Disposable

**Disposable (500 ml / 1000 ml / 3000 ml)**

**Transparent Netting**

- Allows clear view on the fluid bag and the fluid level

**Pocket Design with hanger/hook\***

- Fluid bag can be replaced quickly and easily

**Pressure Gauge**

- Colour coded with overpressure valve
- 360° visible

**Hand Inflator**

- Large volume allows quick inflation

**Three-Way Stop Cock**

- For pressure control
- For quick deflation

**For single use**

- Reduces the potential for cross-contamination
- Latex Free



Product Code With Female Luer Lock	Product Code With Hand Inflator, Stop Cock And Colour Coded Pressure Indicator	Product Code With Hand Inflator And Stop Cock	Size	QTY
56-01-050	56-02-050	56-03-050	500ml	Box 5
56-01-100	56-02-100	56-03-100	1000ml	Box 5
56-01-300	56-02-300	-	3000ml	Each