



Rectal Catheters

Rectal Catheters 40cm

The Rectal Catheter is used to provide access to the human secretion system through rectum to discharge the waste with maximum comfort to the patient. It's also used to relieve flatulence (the presence of excessive gas in the digestive tract) unresponsive to activity or medications

- Atraumatic, soft, 40cm long, frosted & kink resistant PVC tubing
- Sterile, individual peel pouches
- 2 lateral side eyes for efficient drainage
- Proximal end fitted with universal funnel shaped connector for safe connection to waste bag
- Latex Free



Rectal Catheters - Large side eyes / closed tip

Single Wrapped Product Code	QTY	CH/Diam
AERC1816	Each	16ch ø 5.3mm
AERC1818	Each	18ch ø 6.0mm
AERC1820	Each	20ch ø 6.7mm
AERC1822	Each	22ch ø 7.3mm
AERC1824	Each	24ch ø 8.0mm
AERC1826	Each	26ch ø 8.7mm
AERC1828	Each	28ch ø 9.3mm
AERC1830	Each	30ch ø 10.0mm