



Ryles Tubes

Ryles Tubes

Ryles tube used for aspiration of fluids or gases from the stomach, duodenum or intestinal tract through insertion of the nose or mouth. Ryles tube is used to facilitate intestinal decompression.

- Soft, frosted and kink resistance PVC tubing
- Sterile, individual peel pouches
- 105cm in length
- Color coded funnel and connector for easy identification of size and attach luer adaptor
- With X-ray opaque tip, 4 side eyes and closed end
- Latex free

X-Ray line, tip & Luer	QTY	CH/Diam
AERT2306	Each	6ch ø 2.0mm
AERT2308	Each	8ch ø 2.7mm
AERT2310	Each	10ch ø 3.3mm
AERT2312	Each	12ch ø 4.0mm
AERT2314	Each	14ch ø 4.7mm
AERT2316	Each	16ch ø 5.3mm
AERT2318	Each	18ch ø 6.0mm
AERT2320	Each	20ch ø 6.7mm

