



Disposable Self Led Light Proctoscopes

Disposable Self • led • light anoscopies can be used by the specialist both in clinic and in the operating theatre. Made of a transparent plastic material and supplied in individually sealed bags with led light built-in. Led lighting can be performed by pulling out the transparent strip. Guaranteed same light intensity for two hours usage. Compatible, through adapter, to systems of image acquisition (cameras) ALREADY AVAILABLE IN HOSPITALS (for example operating theatres laparoscopic column). ALL ANOSCOPES WITH BUILT-IN LED LIGHT ARE AVAILABLE IN THE STERILE VERSION.

TECHNICAL FEATURES

- Clear transparent plastic proctoscopes
- Range of four size/design to suit varying clinical needs
- Separable anoscope and led lighting for safe disposal
- Flash free
- Guaranteed atraumatic
- Not rechargeable led lighting for safe disposable
- Removable batteries for correct disposal
- Conforms fully with essential requirements of European Council Directive 93/42/EU
- Latex free
- Not Sterile/Sterile



Product Code	Description	QTY
A.4018L	Proctoscope paediatric/diagnostic flute beaked with led lighting built-in Tube length: 43/56mm Tube diameter: 18mm	Box 25
A.4019L	Proctoscope paediatric/diagnostic trunk-end with led lighting built-in Tube length: 54mm Tube diameter: 18mm	Box 25
A.4022L	Surgical/Examination Proctoscope trunk-end with led lighting built-in Tube length: From 65mm to 79mm Tube diameter: 23mm	Box 25
A.4023L	Self LED Light Flute Beaked Proctoscope Disposable Tube Length 65mm Tube Diameter 23mm	Box 25
A.4024L	Self LED Light Trunk End Proctoscope Disposable Tube Length 65mm Tube Diameter 23mm Box 25	Box 25