



Suction Connecting Tubing

Suction Connecting Tubing

A complete system to connect suction sources to suction waste collection systems, suction catheters, yankauers, suction probes and other suction devices. Pennine offer a range of easily connected, interchangeable products in a choice of tubing length and diameter.

- Sterile double wrapped individual peel pouches
- 'F' (Female) connector fits suction machines and male tapered connections for instrumentation.
- 'M' (Male) connector push fits into the 'F' connector and is normally used to connect to the standard funnel fitted on suction catheters. Can also be used to connect lengths of tubing together for longer lengths.
- Application of 'M' connector with 'VC' (vacuum control) vent converts a standard catheter to vacuum control application.
- F = Female
- FF = Female/Female Connectors
- FFM = Female/Female/Male Connectors
- FFM/VC = Female/Female/Male Connectors with Vacuum Control




Product Code	Type	Length	Inside Diam. mm	CH	QTY
CT-4312	F/F	1.8m	5	24CH	Each
CT-4362	F/F	2.7m	5	24CH	Each
CT-4372	F/F	3.5m	5	24CH	Each
CT-4314	F/F/M	1.8m	5	24CH	Each
CT-4364	F/F/M	2.7m	5	24CH	Each
CT-4334	F/F/M	3.0m	5	24CH	Each
CT-4012	F/F	1.8m	6	26CH	Each
CT-4022	F/F	2.0m	6	26CH	Each
CT-4042	F/F	2.7m	6	26CH	Each
CT-4032	F/F	3.0m	6	26CH	Each
CT-4072	F/F	3.5m	6	26CH	Each
CT-4092	F/F	4.0m	6	26CH	Each
CT-4062	F/F	6.0m	6	26CH	Each
CT-4622	F/F	2.0m	8	30CH	Each
CT-4632	F/F	3.0m	8	30CH	Each
CT-4672	F/F	3.5m	8	30CH	Each
CT-4642	F/F	4.0m	8	30CH	Each
CT-4624	F/F/M	2.0m	8	30CH	Each
CT-4634	F/F/M	3.0m	8	30CH	Each
CT-4674	F/F/M	3.5m	8	30CH	Each
CT-4634/VC	F/F/M/VC	3.0m	8	30CH	Each



Suction Connecting Tubing

Flexi-Rib Tubing


- Ribbed connecting tubing for enhanced drapability
- Non-slip outer rib reduces friction
- Resistant to coiling, kinking
- Maximum pressure 500mm
- Multi-fit to all suction devices and sources



Product Code	Type	Length	Inside Diam. mm	CH	QTY
CT-4922	F/F	2.0m	6	26CH	Each
CT-4932	F/F	3.0m	6	26CH	Each
CT-4972	F/F	3.5m	6	26CH	Each
CT-4924	F/F/M	2.0m	6	26CH	Each
CT-4934	F/F/M	3.0m	6	26CH	Each
CT-4722	F/F	2.0m	7	28CH	Each
CT-4732	F/F	3.0m	7	28CH	Each
CT-4772	F/F	3.5m	7	28CH	Each
CT-4724	F/F/M	2.0m	7	28CH	Each
CT-4734	F/F/M	3.0m	7	28CH	Each
CT-4774	F/F/M	3.5m	7	28CH	Each
CT-4714/VC	F/F/M/VC	1.8m	7	28CH	Each
CT-4741/VC	F/F/M/VC	2.5m	7	28CH	Each

Link Suction Connecting Tubing/Bubble Tubing

The Pennine Link system is used for connecting suction sources to suction waste collection systems, suction catheters, yankauer suckers, suction probes and other suction devices. Tubing extruded with an integral female funnel and a tapered male end to provide adaptable connections. The tube can be modified by cutting with scissors.



Product Code	Type	Length	Inside Diam. mm	CH	QTY
CT-4513	Plain F/F	1.8m	5	24CH	Each
CT-4543	Plain F/F	2.7m	5	24CH	Each
CT-4511	Plain M/F	1.8m	5	24CH	Each
CT-4541	Plain M/F	2.7m	5	24CH	Each
CT-4713	Plain F/F	1.8m	7	28CH	Each
CT-4783	Plain F/F	3.7m	7	28CH	Each
CT-4741	Plain M/F	2.7m	7	28CH	Each
CT-4743	Plain F/F	2.7m	7	28CH	Each
CT-4711	Plain M/F	1.8m	7	28CH	Each