










Tourniquet Cuffs Accessories

Tourniquet Sleeves

- Protects soft tissues underneath the cuff
- Reduces wrinkling, pinching and shearing of soft tissue
- Sterile and non sterile version
- For single use

Patient	Size	Color code	Product Code Clean packed	Product Code Sterile packed	Qty
Infant	≤ 20 cm		23-00-700	23-00-700-1	Box 10
Child	≤ 30 cm		23-00-710	23-00-710-1	Box 10
Arm	≤ 35 cm		23-00-711	23-00-711-1	Box 10
Arm long	≤ 46 cm		23-00-712	23-00-712-1	Box 10
Leg	≤ 61 cm		23-00-722	23-00-722-1	Box 10
Leg long	≤ 76 cm		23-00-727	23-00-727-1	Box 10
Leg extra long	≤ 86 cm		23-00-728	23-00-728-1	Box 10




Cuff Cover

Keeps the cuff free from soiling and allows sealing the cuff with adhesive tape without damaging the cuff

- For Single Use


Size	Product Code Clean packed	Product Code Sterile packed	Qty
Small	20-11-111	20-11-111-1	Box 50
Medium	20-11-333	20-11-333-1	Box 50
Large	20-11-222	20-11-222-1	Box 50



Measuring Tape

Colour coded, for easy identification of patient extremity and optimal cuff selection

Version	Product Code	Qty
Single use	20-95-222	Box 50





Roll-On Cuffs

Allow fast and easy blood exsanguination within seconds (an alternative to Esmarch Bandage) autoclavable, Material: Silicone, Colour: transparent



Product Code	Description	QTY
21-95-400	Roll-On Cuff Set, small complete in case Consisting of: • Roll-On Cuffs # 1-4 • Manometer with connecting tube • Measuring tape • Brake-block, large or small	Each
21-95-401	Roll-On Cuff Set, large complete in case consisting of: • Roll-On Cuffs # 3-6 • Hand pump with manometer and connecting tube • Measuring tape • Brake-block, large or small	Each



Product Code	Description	QTY
21-95-712	Roll-On Cuffs Size2 26-33cm Light blue	Each
21-95-719	Roll-On Cuffs Size3 31-40cm Yellow	Each
21-95-722	Roll-On Cuffs Size4 38-49cm Red	Each
21-95-727	Roll-On Cuffs Size5 47-60cm Green	Each
21-95-729	Roll-On Cuffs Size6 58-70cm Dark blue	Each



Product Code	Description	QTY
21-95-350	Hand Pump with manometer and connecting tube (100 cm)	Each
21-95-300	Manometer with connecting tube (100 cm)	Each
21-95-302	Silicone connection tubing 100 cm long, with positive locking connector and inflation needle	Each
21-95-200	Brake-Block reusable, Material: Silicone, Small	Each
21-95-201	Brake-Block reusable, Material: Silicone, Large	Each
21-95-111	Measuring Tape color coded, for easy identification of patient extremity and cuff size	Box 10



BOSCO
MEDICAL AUSTRALIA

VBM Medizintechnik

Esmarch Bandages

- Allows safe exsanguination of any extremity
- Large range of sizes available in 5 different widths
- Reusable version
- Single use version, sterile



REUSABLE

Material: Polyisoprene (PI)

**LATEX
FREE**

Product Code	Description	QTY
20-20-060	6cm x 5m Length	Each
20-20-080	8cm x 5m Length	Each
20-20-100	10cm x 5m Length	Each
20-20-150	15cm x 5m Length	Each



SINGLE USE

Material: Polyisoprene (PI)

**LATEX
FREE**

Product Code	Description	QTY
20-21-060-1	6cm x 4m Length	Box 10
20-21-080-1	8cm x 4m Length	Box 10
20-21-100-1	10cm x 4m Length	Box 10
20-21-150-1	15cm x 4m Length	Box 10



BOSCO
MEDICAL AUSTRALIA

Bosco Medical Australia Pty Ltd
Unit 23, 160 Lytton Road
Morningside, QLD 4170

P: + 61 7 3399 1799 F: + 61 7 3395 6736
Email: orders@boscomed.com.au
Website: www.boscomed.com.au