



Humidifiers & Water Traps

For dependable, cost effective performance, Salter humidifiers deliver a durable, disposable product that's ideal for use with all oxygen sources; concentrators, liquid tanks, or wall outlets. Unique 360° rectangular diffusion ports offer quiet operation with improved hydration and less system back pressure.



Standard Bubble Humidifiers

Available with either 3 PSI or 6 PSI pop off alarms

Description	QTY	Product Code
Bubble Humidifier 350ml, 3psi pop off (brown lid)	Box 50	7100-0-50
Bubble Humidifier 350ml, 6psi pop off (black lid)	Box 50	7600-0-50

High Flow Bubble Humidifiers

Ideal for use with high-output oxygen concentrators, wall outlet source gas, or other applications where flow rates of 6 to 15 LPM are required and an audible pressure relief alter is desired. Effectively pairs with Salter 1600HF and 16SOFT-HF nasal cannulas.

Description	Product Code	QTY
Bubble Humidifier 350ml, (2-tone, beige and black lid)	7900-0-25	Box 25

Salter Inline Water Traps

A cost-effective solution for your long term patients: simple, durable construction, with design features that minimize the hazard of snagging on furniture or other obstructions.



Description	QTY	Product Code
Water trap, clear body with self-sealing cap	Box 25	7000-0-25
Water trap, clear body with threaded end cap	Box 25	7001-0-25