



SnapIT Ampoule Opener

Clicksmart SnapIT is a simply safer way to open glass ampoules. This easy-to-use multi-use ampoule opener avoids the unacceptable painful and dangerous sharps injury rates seen when ampoules are opened by hand.

SnapIT provides an essential safety solution to:

- Provide best practice occupational health and safety
- Avoid infection and contamination risk caused by sharps injuries
- Reduce time lost to sharps injuries
- Ensure a safe and effortless ampoule opening every time due to the extended lever action
- Prevent wastage of expensive medication contained in ampoules

AVOID THIS TYPE OF INJURY




Snap it, Simply and Safely


With SnapIT you can safely:

1. Insert the ampoule lid into your multi-use SnapIT
2. Snap open the ampoule
3. Eject the ampoule lid directly into the sharps container

SnapIT avoids ampoule sharps injuries by keeping both hands away from the sharp edges during ampoule opening. One hand holds the base of the ampoule while the other hand simply uses the SnapIT to snap off the ampoule lid. The ampoule lid is contained safely within the SnapIT until ejected.



Causative Events for Needlestick and Sharps Injuries³



Product Code	Description	QTY
SN-01RB	SnapIT Ampoule Breaker, Personal, Blue	Each
SN-01RP	SnapIT Ampoule Breaker, Personal, Purple	Each
SN-02L	SnapIT Ampoule Breaker, Lite, Large	Each
SN-02R	SnapIT Ampoule Breaker, Lite, Regular	Each
TE-01L	SnapIT Ampoule Breaker, Trolley Version, Large	Each
TE-01R	SnapIT Ampoule Breaker, Trolley Version, Regular	Each
TE-01XL	SnapIT Ampoule Breaker, Trolley Version, XLarge, Blue	Each